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Enhancing Class II treatment efficiency with Clarity™ SL Self-Ligating Brackets and Forsus™ Class II Correctors

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Introduction

The Class II malocclusion is considered one of the most challenging treatment needs to orthodontists. Besides the skeletal and/or dental characteristics, facial features can also influence treatment options. In this case report, a half-cuspid Class II malocclusion patient was successfully treated with Forsus™ Class II Correctors, enabling reduction in treatment time and patient's compliance. Clarity™ SL Self-Ligating Brackets were selected to improve aesthetic appearance and performance (low friction during initial phases and control during Class II correction). The MBT™ System prescription was essential to achieve excellent torque control, especially in the lower anteriors.

Patient

Female; 12 years, 11 months

Patient's Main Concern

Misalignment of anterior teeth and spacing between upper central incisors.

X-ray Findings

- Complete permanent dentition
- Roots mineralization beginning on third molars
- Short roots noted on teeth #11 and 21

Dental Analysis

- Half-cusp Class II Division 2 malocclusion
- Deep overbite (5 mm)
- Upper central incisors extruded and tipped distally with a 2 mm diastema
- Upper lateral incisors tipped labially
- Lower posterior teeth tipped lingually
- Accentuated lower curve of Spee (4 mm)

Treatment Plan

- Upper/Lower Clarity SL ceramic self-ligating brackets -.022" slot – MBT System Prescription
- Molar tubes bonded to upper and lower molars -.022" slot – MBT Prescription

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- Buildups on lower first molars tubes to open the bite allowing lower brackets positioning
- Alignment and leveling of upper and lower arches
- Forsus™ Appliance for Class II correction (bonded first upper molar tubes were replaced with bands)

Treatment		20 months (February 2013 – October 2014)	
Mx	February 2013	Indirect	.014" SE.016" SE.014" and .016" SE.017" x.025" SE.019" x.025" SE.019" x.025" SS.019" x.025" Braided
Md	April 2013	Indirect	.014" SE.016" SE.018" SE, .017" x.025" SE.019" x.025" SE.019" x.025" SS.019" x.025" Braided
Forsus™ Treatment	4 months (July 2013 – November 2013)		
# of visits	18		
Emergencies	1		

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Retention

- Removable appliance with bite plane on upper arch to avoid lower incisors extrusion
- 3×3 fixed lingual retainer on lower arch

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




Doctor's Notes

1. Lower brackets and tubes were bonded in the third month of treatment after an initial alignment of upper incisors. Buildups in first lower molar were planned to prevent upper teeth enamel wear and lower brackets failure.
2. First lower molars buildups were removed progressively during treatment allowing extrusion of these teeth during lower curve of Spee leveling.
3. SE denotes NiTi Super Elastic wire and SS denotes stainless steel wire.
4. This case was treated with Forsus Correctors, adapted to bands on first upper molars. Currently, it is recommended to use the Forsus™ Wire Mount that requires no bands.
5. Forsus Correctors were adapted to lower canine's brackets objecting a more horizontal vector during Class II correction due patient's horizontal skeletal pattern. Adapting Forsus Correctors on first premolar's bracket could produce a more vertical vector tending to intrude upper first molar that was not a desirable effect in this case.
6. The effects of Class II treatment with Forsus Correctors were in accordance to the reported in previous studies, including distal movement of upper molars, mesial movement of lower molars and proclination of lower incisor. No skeletal effect was observed.
7. After Forsus Appliance correction, nocturnal Class II elastic were applied (1/4 in./8 oz.). During finishing some brackets were repositioned to improve occlusal relationship.
8. MBT System incisors torque prescription (U1=17°; U2=10°; L1 and L2=-6°) allowed an excellent torque control of these teeth especially in the lower arch. Proclination of lower incisors was well controlled despite the effects of the Forsus Corrector therapy and Curve of Spee leveling.

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Some article content was previously published in „The Tandem Archwire Concept with Self-Ligating Brackets.“, Journal of Clinical Orthodontics, 2014; Vol XLVIII, No 4. Used with permission.
3M, Clarity, Forsus, and SmartClip are trademarks of 3M. Used under license in Canada.
© 3M 2016. All rights reserved.